

For Trial Patients


How to Use Your HFX App

How to download the HFX App

1. Open your mobile device from your Home screen. Tap the blue icon that says **App Store**.
2. Tap the **Search** icon on the bottom right side of your screen.
3. In the search bar type, “**HFX iQ**” and then tap **Search**.
4. Tap on the **HFX App** and tap **GET** to download the app.

Note

You will need to know your existing Apple ID and password to download the HFX App from the Apple App Store.



Or scan this code to go to the Apple App Store.

Get familiar with your Home screen

1. Stim ON/OFF Button
2. Therapy Settings
3. Change Program Button
4. Start Assessment
5. Battery Indicator
6. Connection Status
7. Menu Tab



For Trial Patients

How to Use Your HFX App

How to change programs

1. Tap the **Change Program Button** on the lower right of the circle.
2. To change your program setting, scroll through the program carousel until the circle is over the program you want.
3. To change your level setting, scroll through each level on the level wheel until the level you want is selected.
4. Tap **Done** to save your changes.




Trial guidelines

- Keep your mobile phone within close distance to your trial device because the two will communicate.
- You will make all changes using your HFX App.
- Check the Stim ON/OFF Button on your HFX App twice a day to confirm your trial device is ON.

How to complete an assessment

- When you're within your assessment window, the Start Assessment button will turn yellow.



Your Nevro HFX Care Team Contact Information

1-844-331-1001

Monday - Friday 9:00am - 5:00pm



Nevro Corp.
1800 Bridge Parkway Redwood City, CA 94065
1.650.251.0005 • 1.888.956.3876

HFX, the HFX logo, HFX iQ, the HFX iQ logo, NEVRO, and the NEVRO logo are trademarks or registered trademarks of Nevro Corp.
© 2023 Nevro Corp. All Rights Reserved.
S001632 V1.0
ML-00035 Rev. A